

MONOSCREEN[®] Ab ELISA

Instructions for use
 BIO K 466-Besnoitia besnoiti_NO_(EN)_V02
 15/03/2022

Monoscreen AbELISA *Besnoitia besnoiti*

Reference : BIO K 466

ELISA Kit for serodiagnosis of Besnoitiosis

Monowell, blocking

In vitro and strictly veterinary use



| Sample | Species | Individual analysis | Pool analysis*, possible up to |
|-----------------------|---------|---------------------|--------------------------------|
| Blood sera and plasma | Cattle | ✓ | 10 |

* This is done in accordance with the legislation in force in your country, the certifying body or the recommendations made by the NRL when they exist. Mixtures must be made volume to volume, i.e. by taking the same volume of each of the sera making up the mixture.


Presentation

| Product reference | BIO K 466/2 |
|-------------------|----------------------------|
| Format | 2 plates, strip of 8 wells |
| Reactions | 192 tests |

Composition of the kit

| Provided material | BIO K 466/2 |
|--------------------------------|-------------|
| Microplates | 2 |
| Washing solution (20X) | 1 X 100 mL |
| Colored dilution solution (1X) | 1 X 60 mL |
| Conjugate (50X) | 1 X 0,55 mL |
| Positive control | 1 X 0,5 mL |
| Negative control | 1 X 0,5 mL |
| Single component TMB (1X) | 1 X 25 mL |
| Stop solution (1X) | 1 X 15 mL |

Revision history

| Date | Version | Modifications |
|--|---------|--|
| 07/03/2022 | V01 | First version |
| 15/03/2022  | V02 | Addition of columns "individual analysis" and "pool analysis*, possible up to" |

Note : minor changes to typography, grammar and formatting are not included in the revision history

A. Introduction

Besnoitia besnoiti, the causative agent of bovine besnoitiosis, is an obligate intracellular protozoan. The disease affects mainly young cattle. Besnoitiosis is epizootic in the south of France, but is now widely distributed in Africa, Asia and in Southwestern Europe. The most likely pathway of transmission would be transcutaneous, by stinging insects (tabanids, stomox).


During infection, an incubation stage of 3 to 6 days is followed by 3 successive clinical stages:

- A febrile stage of 3 to 7 days; the tachyzoites multiplication in endothelial cells of blood vessels increases the animal temperature
- A second stage of 1 to 2 weeks; the bradyzoites cysts generate subcutaneous oedema
- A chronic stage of several months, characterized by alopecia and scleroderma. The skin becomes then markedly thickened and wrinkled, and parasitic cysts are observed on conjunctiva and sclera. This ultimate phase leads generally to the death of the animal or to its euthanasia.

Serologic tests are available for the detection of the specific antibodies of *Besnoitia besnoiti* present in the chronic stage. In order to avoid the transfers contaminated animals and to control the spread of the bovine besnoitiosis, it is essential to use diagnostic tools that detect the pathogens at the early stages of the disease.

B. Test principle

The test uses 96-well microtitration plates sensitized by *Besnoitia besnoiti*'s culture lysate. The operator deposits the previously diluted test sera in the microplate's wells. After 45 minutes' incubation and a washing step, the operator adds the conjugate, which is a specific monoclonal antibody against *Besnoitia besnoiti* coupled to peroxidase. After incubating and washing the preparation, the operator adds the chromogen tetramethylbenzidine (TMB). This chromogen has the advantage of being more sensitive than the other peroxidase chromogens and not being carcinogenic. The intensity of the color is inversely proportionate to the sample's serum titer. Positive and negative control sera are provided with the kit to be able to validate the test results.



C. Additional material and required equipment (not provided)

- Distilled/demineralized water
- Graduated mono- or multichannel pipettes (2-20 μ L, 20-200 μ L et 100-1000 μ L range) and single-use tips
- Microplate reader (450nm filter)
- Microplate washer
- Incubator at 37 \pm 2 $^{\circ}$ C
- Standard laboratory equipment: graduated cylinder, tube rack, lid, ...

D. Precautions for use

- The reagents must be kept between +2 et +8 $^{\circ}$ C.
- Unused strips must be stored with the desiccant in the hermetically sealed aluminum envelope.
- Do not use reagents beyond shelf-life date.
- Do not use reagents from other kits.
- Make sure to use distilled/demineralized water.

- The stopping solution contains 1 M phosphoric acid. Handle it carefully.
- Used material must be disposed of in compliance with the legislation in force regarding environmental protection and biological waste management.
- Keep the TMB solution away from light.

E. Preparation of solutions

- The solutions are to be prepared extemporaneously.
- The washing solution must be diluted 20-fold in distilled/demineralized water. The cold solution crystallizes spontaneously. Bring the vial to 21 \pm 3 $^{\circ}$ C to make sure that all crystals have disappeared; mix the solution well and withdraw the necessary volume.
- The dilution solution is ready to use. The dilution solution is colored in yellow. It is used for dilution of samples, positive and negative serums, and conjugate.
- The conjugate must be diluted 50-fold in the dilution solution.
- The stopping solution is ready to use.
- The TMB solution is ready to use. It must be perfectly colorless.

F. Procedure

- Bring all the reagents to 21 \pm 3 $^{\circ}$ C before use.
- Carefully read through the previous points.

1. Distribute 50 μ L/well of dilution solution 1X. Add 50 μ L/well of serum samples and of positive and negative controls. Homogenize by pipetting up and down.

N.B.: To avoid differences in incubation time between samples, sample dilutions and reference dilutions can be prepared in a dilution microplate before transfer (200 μ L) into the test microplate using a multi-channel pipette.

2. Cover and incubate the plate at 37 \pm 2 $^{\circ}$ C during 120 \pm 5 min.
3. Remove the content of the microplate. **Wash the microplate 3 times with 300 μ L of washing solution per well.** Avoid the formation of bubbles in the wells and the desiccation of the microplate between each wash.
4. Add **100 μ L of diluted conjugate** per well. Cover with a lid and incubate the plate at 37 \pm 2 $^{\circ}$ C during 30 \pm 2 min.
5. Remove the content of the microplate. **Wash the microplate 3 times with 300 μ L of washing solution per well.** Avoid the formation of bubbles in the wells and the desiccation of the microplate between each wash.
6. Distribute **100 μ L of TMB solution** per well.
7. Incubate at 21 \pm 3 $^{\circ}$ C during 10 \pm 1 min away from the light, without covering.
8. Distribute the stopping solution at rate of **50 μ L per well.** Color changes from blue to yellow.
9. Record the optical densities using a plate spectrophotometer with a 450 nm filter within 5 minutes after adding the stopping solution.

G. Validation of results

The test can only be validated if :

- the difference between positive and negative serum optical density readings is greater than 0,700.

$$OD_{\text{negative serum}} - OD_{\text{positive serum}} > 0,700$$

- the positive serum's inhibition percentage (%inh) is greater than 50%.

$$\%inh_{\text{positive serum}} > 50\%$$

H. Interpretation of results

Calculate for each sample its inhibition percentage (%inh) using the following formula :

$$\%inh = \frac{OD_{\text{negative serum}} - OD_{\text{sample}}}{OD_{\text{negative serum}}} * 100$$

| | Results | Status |
|--------|-------------|----------|
| Sample | %inh < 40% | Negative |
| | %inh ≥ 40 % | Positive |

Get the interpretation of your results quickly and easily using **AnalysiScreen**, our free online platform, available on our website : <https://www.biox.com>





AnalysiScreen™ is the new module for reading and interpreting all types of Monoscreen™ and Multiscreen™ ELISA plates. AnalysiScreen™ is:

- Free
- Accessible online via our website: <https://www.biox.com>
- Updated in real time
- Compatible with all Bio-X Diagnostics plate designs
- Very easy to use



SCAN ME

50 μ L of dilution solution 1X + 50 μ L of serum
 50 μ L of dilution solution 1X + 50 μ L of positive serum
 50 μ L of dilution solution 1X + 50 μ L of negative serum

Add 100 μ L of diluted conjugate (dilution 1/50)



Add 100 μ L of TMB



Add 50 μ L of stopping solution

Record the optical densities



* The notes do not replace the instructions for use of which they are a synthesis.



Distribuito in ITALIA da
Li StarFish S.r.l.
 Via Cavour, 35
 20063 Cernusco S/N (MI)
 telefono 02-92150794
 info@listarfish.it
 www.listarfish.it