



Distribuito in ITALIA da  
**Li StarFish S.r.l.**  
Via Cavour, 35  
20063 Cernusco S/N (MI)  
telefono 02-92150794  
fax 02-92157285  
info@listarfish.it  
www.listarfish.it

GLOBAL DIAGNOSTICS B



# MICRO SHAKE 1200


## ELISA Plate Shaker




**Synchronised Harmony with Skaking & Incubation**

## FEATURES


- Simultaneous orbital shaking & incubation operation
- Buzzer indication on completion of incubation
- Indication of Remaining Time
- Display of temperature on pressing TEMP Key



RPM  
400 to 700



Incubator



Shaker

## SPECIFICATIONS

Operating Modes	:	Shaker Incubator Shaker & Incubator
Temperature Control		
Temperature Range	:	37° C to 42° C
Resolution	:	1° C
Incubation Time	:	1 to 999 min
Shaker		
Frequency	:	400 to 700 RPM
Operating Position	:	On horizontal flat, rigid and vibration free surface
Operating Conditions		
Temperature	:	From + 18° C to 32° C
Relative Humidity	:	Up to 80%
Storage Conditions		
Temperature	:	From - 10° C to 35° C
Relative Humidity	:	Up to 80%
Enclosure	:	ABS Fire retardant
Size (cm)	:	28 X 25 X 16 (l X b X h)
Weight	:	3 Kg (Approx)

*Due to continuous development, specifications are subject to change without prior notice.*



**Global Diagnostics B**

Vlasmeeer 5, 2440 Mol, BELGIUM,  
Tel.: +32468220039  
info@gdms.be