



AutoMag Solution



Rapid implementation with prevalidated and preprogrammed script.

Optimized extraction kits, plastic consumable = **complete** solution.

Superparamagnetic beads produced on-site and used in all of our extraction kits.

Easy to use prefilled kits and UV lamp prevents the risk of cross-contamination, improving user safety and product quality.

We able to custom your protocol

- Easy implementation of a new protocol by USB stick memory
- Easy to install and to move without external intervention

Solid and without liquid metrology : the particles are moved to the different wells using the magnet head and avoid the contamination by aerosols

Li StarFish distribuisce:



PRODOTTI PER LA RICERCA MEDICA

- Compact : size of a 50cm cube
- Scripts are provided updates to run for different applications
- Free updates

Li StarFish distribuisce:

