



Human Parathyroid Hormone Related Protein ELISA Kit

This sandwich kit is for the accurate quantitative detection of human Parathyroid Hormone Related Protein (also known as PTHrP) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Cat.No BTB-E1018Hu

Standard Curve Range: 8ng/L – 1500ng/L

Sensitivity: 4.59ng/L

Storage: The kit could be stored for 6 months under the temperature of 2-8°C and could be stored for longer time under the temperature below -20°C. The ELISA plate is packed in vacuum package and the strip is dismountable. If it's not used up once after the dismounting, put it in the provided seal bag and put in the drying agent and store it under the temperature of 2-8°C. It could be used again within one month. Prevent the repeated frost-thawing for the reagent.

Kit Contents

- Standard Solution
- Pre-coated ELISA Plate
- Standard Diluent
- Streptavidin-HRP
- Stop Solution
- Substrate Solution A
- Substrate Solution B
- Wash Buffer Concentrate (25x)
- Biotin-Conjugate Antibody
- User Instruction
- Plate Sealer
- Zipper bag

Li StarFish distribuisce: