



PRODOTTI PER LA RICERCA MEDICA

Li StarFish S.r.l.
Via Cavour, 35 - 20063 Cernusco S/N (MI), Italy
Tel. +39-02-92150794 - Fax. +39-02-92157285
info@listarfish.it - www.listarfish.it

Monoclonal antibody to Human Interleukin-10

Antibody:	9D7
Product code:	3430-3-1000
Quantity:	1000 µg
Immunogen:	Recombinant human IL-10
Isotype:	Rat IgG1
Specificity:	Native and recombinant human IL-10
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	For quantification of human IL-10 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-10 producing cells using ELISpot. 9D7 is recommended as coating mAb in combination with detection mAb 12G8 (product code 3430-6). Immunostaining of human IL-10.
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

Note; for research use only.

Li StarFish distribuisce:

