

Monoclonal antibody to Mouse TNF- α

BIOTINYLATED

Antibody:	MT11B10
Product code:	3511-6-250
Quantity:	250 μ g
Immunogen:	Recombinant mouse TNF- α
Isotype:	Rat IgG2a
Specificity:	Native and recombinant mouse TNF- α
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	For quantification of TNF- α in solution e.g. cell culture supernatants using ELISA and for enumeration of TNF- α producing cells using ELISpot. MT11B10 is recommended as detection mAb in combination with coating mAbs MT1C8/23C9 (product code 3511-3).



Distribuito in ITALIA da
Li StarFish S.r.l.
Via Cavour, 35
20063 Cernusco S/N (MI)
telefono 02-92150794
fax 02-92157285
info@listarfish.it
www.listarfish.it

Note; for research use only.