

Supplier	Ademtech
Instrument	AutoMag Solution PRIME
Technology	<ul style="list-style-type: none"> • Magnetic beads • Magnetic Rods • Plate to fill and/or prefilled reagents plate • Intuitive software flexible parameters (mixing speed, temperature, time) • Protocol development achievable by Ademtech
Number of samples processed	<ul style="list-style-type: none"> • 12 samples in parallel • Up to 24 samples sequentially (2 plates)
Number of steps possible per plate	<ul style="list-style-type: none"> • 8
Consumables	<ul style="list-style-type: none"> • Microtiter Deepwell 96 plate: 50 pc • 12- tip comb for 96 DW plate: 50 pc • Elution strip: 40 pc • Cap for elution strip: 40 pc
Price of Consumables 2018 (End User)	<ul style="list-style-type: none"> • Microtiter Deepwell 96 plate : 202€ • 12- tip comb for 96 DW plate : 67€ • Elution strip : 53€ • Cap for elution strip : 40€
Applications	<ul style="list-style-type: none"> • Sample preparation • MeDIP • HMeDIP • ChIP • Re-ChIP • DNA purification (ChIP DNA Prep)
Size, dimension	<ul style="list-style-type: none"> • 40 X 46 X 35 cm • 16 kg
Process time	<ul style="list-style-type: none"> • ChIP protocol : mini 1h30 (ChIP Flash)
Working volume	<ul style="list-style-type: none"> • 50 µl-1000 µl (96 DW plate, 12 pin magnet head) • 30-130 µl (elution strip, 12 pin magnet head)
Heating system	<ul style="list-style-type: none"> • Block temperature from +10°C to +75°C in RT
Heating temperature for elution strip	<ul style="list-style-type: none"> • Block temperature from +4°C to +75°C in RT



Distribuito in ITALIA da
Li StarFish S.r.l.
 Via Cavour, 35
 20063 Cernusco S/N (MI)
 telefono 02-92150794
 fax 02-92157285
 info@listarfish.it
 www.listarfish.it