



Li StarFish S.r.l.
Via Cavour, 35
20063 Cernusco S/N (MI), Italy
Tel. +39-02-92150794 - Fax. +39-02-92157285
info@listarfish.it - www.listarfish.it

Anti-MMAF Polyclonal Antibody

Catalog no.	EPI-30700
Target	Monomethyl auristatin F
Host Species	Mouse
Isotype	IgG
Clonality	Polyclonal
Purification Method	Affinity Chromatography
Conjugation	Unconjugated
Concentration	25 µg/mL
Validated Application	ELISA
Buffer	0.01M PBS, 0.1% BSA, Proclin-300
Storage Temp	-20 °C (-80 °C for long term)
Note:	This antibody comes as 50x concentrated and is used in the MMAF-ADC competitive EIA (Cat. no. EPI-KTR-758)

Li StarFish distribuisce:

