



Li StarFish S.r.l.
Via Cavour, 35
20063 Cernusco S/N (MI), Italy
Tel. +39-02-92150794 - Fax. +39-02-92157285
info@listarfish.it - www.listarfish.it

Anti-MMAE Antibody

Catalog no.:	EPI-MAB30699
Target:	Monomethyl auristatin E
Host Species:	Mouse
Isotype:	IgG 1
Clonality:	Monoclonal
Clone no.:	2E2
Purification Method:	Affinity Chromatography
Conjugation:	Unconjugated
Concentration:	0.5 - 2.0 mg/mL
Validated Application:	ELISA
Buffer:	0.01M PBS
Storage Temp:	-20 °C (-80 °C for long term)

For research use only.

Not for use in diagnostics procedures.

Li StarFish distribuisce:

