



Li StarFish S.r.l.
Via Cavour, 35
20063 Cernusco S/N (MI), Italy
Tel. +39-02-92150794 - Fax. +39-02-92157285
info@listarfish.it - www.listarfish.it

Anti-MMAF Antibody

Catalog no.:

EPI-MAB30700

Target:

Monomethyl Auristatin F

Host Species:

Mouse

Isotype:

IgG 1

Clonality:

Monoclonal

Clone no.:

3D6

Purification Method:

Affinity Chromatography

Conjugation:

Unconjugated

Concentration:

0.5 - 2.0 mg/mL

Validated Application:

ELISA

Buffer:

0.01M PBS

Storage Temp:

-20 °C (-80 °C for long term)

For research use only

Not for use in diagnostics procedures.

Li StarFish distribuisce:

