



## Recombinant Human PTHrP

**Catalogue Number:** PEP-100-09

**Description:**

PTHrP is a polypeptide hormone produced by almost every tissue of the body. PTHrP is closely related to parathyroid hormone (PTH), which is secreted from the parathyroid gland, and plays a central role in regulating the extracellular concentrations of calcium and phosphorous. Recombinant Human PTHrP is a 9.8 kDa linear polypeptide of 86 amino acid residues.

**Source:** *E.coli*

**Synonyms:** Parathyroid Hormone-related Protein

**AA Sequence:** AVSEHQLLHD KGKSIQDLRR RFFLHHLIAE IHTAEIRATS EVSPNSKPSP NTKNHPVRFG SDDEGRYLTQ ETNKVETYKE QPLKTP

**Purity:** ≥ 98% by SDS-PAGE gel and HPLC analyses.

**Biological Activity:** Not available.

**Cross Reactivity:**

Human , Mouse

**Country Of Origin:** USA

Not for human use.

*Li StarFish distribuisce:*